



Product Datasheet

BTNL8 antibody (orb35306)

Description	Rabbit polyclonal antibody to BTNL8. The protein is a single-pass type I membrane protein that belongs...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	BTNL8
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56748
Uniprot ID	Q6UX41
NCBI	NP_001035552.1 , NP_001153181.1 , NP_079126.1 , NP_001153179.1 , NP_001153182.1 , NP_001153180.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of in MDA-MB435 cells...