

Product Datasheet

NCS1 antibody (orb353077)

Description	NCS1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Neuronal calcium sensor 1 protein (1-190AA)
Target	NCS1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62166
Dilution Range	WB:1:500-1:2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.