



Product Datasheet

FATE1 antibody (orb35308)


Description	FATE1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Fetal and adult testis-expressed transcript protein (1-162AA)
Target	FATE1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q969F0
Dilution Range	IHC-P:1:100-500, IF:1:50-500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...



Immunohistochemical staining of human ov...



Immunohistochemical staining of human te...