


Product Datasheet

ALDH1A2 antibody (orb353081)

Description	ALDH1A2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Retinal dehydrogenase 2 protein (1-180AA)
Target	ALDH1A2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O94788
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...