



Product Datasheet

CREB3L2 antibody (orb353121)

Description	CREB3L2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Cyclic AMP-responsive element-binding protein 3-like protein 2 protein (1-247AA)
Target	CREB3L2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q70SY1
Dilution Range	WB:1:200-1000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis
of mouse liver & l...



Immunohistochemical
staining of human
ki...