

Product Datasheet

AMBRA1 antibody (orb353144)

Description	AMBRA1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Activating molecule in BECN1-regulated autophagy protein 1 protein (1-240AA)
Target	AMBRA1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9C0C7
Expiration Date	12 months from date of receipt.