



Product Datasheet

MAP2K3 antibody (orb353154)


Description	MAP2K3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Dual specificity mitogen-activated protein kinase kinase 3 protein (1-170AA)
Target	MAP2K3
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P46734
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pl...



Immunohistochemical staining of human to...



Immunohistochemical staining of human en...