




Product Datasheet

CMTM7 antibody (orb35317)

Description	Rabbit polyclonal antibody to CMTM7. The protein is also referred as CKLFSF7, and belongs to the che...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	CMTM7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18834
Uniprot ID	Q96FZ5
NCBI	NP_852137.1 , NP_612419.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of in Ramos cell l...



Western blot analysis of 293 cell lysate...