

Product Datasheet

CDK16 antibody (orb353214)

Description	CDK16 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Cyclin-dependent kinase 16 protein (277-496AA)
Target	CDK16
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q00536
Dilution Range	WB:1:200-1000
Expiration Date	12 months from date of receipt.

