


Product Datasheet

STXBP4 antibody (orb353255)

Description	STXBP4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Syntaxin-binding protein 4 protein (344-553AA)
Target	STXBP4
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6ZWJ1
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of PC-3 whole cell...