


Product Datasheet

CAMK4 antibody (orb353290)

Description	CAMK4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Recombinant Human Calcium/calmodulin-dependent protein kinase type IV protein (304-473AA)
Target	CAMK4
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16566
Dilution Range	WB:1:1000-5000, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunoprecipitation
analysis of Jurkat
W...



Western blot
analysis of 293T
whole cel...