


Product Datasheet

CD200R1 antibody (orb353291)

Description	CD200R1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Cell surface glycoprotein CD200 receptor 1 protein (226-325AA)
Target	CD200R1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TD46
Dilution Range	WB:1:200-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse lung tiss...