



Product Datasheet

GNB1 antibody (orb35330)

Description	Rabbit polyclonal antibody to GNB1. The protein is a member of the WD repeat G protein beta family, ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	GNB1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37377
Uniprot ID	P62873
NCBI	NP_002065.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of in Jurkat cell ...