



Product Datasheet

IL7R antibody (orb353379)


Description	IL7R antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Interleukin-7 receptor subunit alpha protein (21-239AA)
Target	IL7R
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16871
Dilution Range	WB:1:500-1000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of K562(lane 1), J...



Western blot analysis of hepG2 cell usin...



Immunohistochemical staining of human co...