

Product Datasheet

CHRM3 antibody (orb353382)

Description	CHRM3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Muscarinic acetylcholine receptor M3 protein (270-470AA)
Target	CHRM3
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20309
Expiration Date	12 months from date of receipt.