



Product Datasheet

GATA3 antibody (orb353407)


Description	GATA3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Trans-acting T-cell-specific transcription factor GATA-3 protein (1-260AA)
Target	GATA3
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P23771
Dilution Range	WB:1:500-2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mcf-7 whole cell...



Immunohistochemical staining of human to...



Immunohistochemical staining of human ep...