



## Product Datasheet

### AKR1C3 antibody (orb353409)


|                            |                                                                            |
|----------------------------|----------------------------------------------------------------------------|
| <b>Description</b>         | AKR1C3 antibody                                                            |
| <b>Species/Host</b>        | Rabbit                                                                     |
| <b>Reactivity</b>          | Human                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                               |
| <b>Tested Applications</b> | ELISA, IHC, IP, WB                                                         |
| <b>Immunogen</b>           | Recombinant Human Aldo-keto reductase family 1 member C3 protein (1-323AA) |
| <b>Target</b>              | AKR1C3                                                                     |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                          |
| <b>Form/Appearance</b>     | Liquid                                                                     |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.              |
| <b>Note</b>                | For research use only                                                      |
| <b>Isotype</b>             | IgG                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                 |
| <b>Uniprot ID</b>          | <b>P42330</b>                                                              |
| <b>Dilution Range</b>      | WB:1:500-5000, IHC-P:1:20-1:200                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                            |



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of K562 whole cel...



Immunoprecipitation analysis of HepG2 wh...