



Product Datasheet

IRF4 antibody (orb353438)


Description	IRF4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Interferon regulatory factor 4 protein (180-451AA)
Target	IRF4
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15306
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human colon tissue showing IRF4 expression.



Immunohistochemical staining of human colon tissue showing IRF4 expression.



Immunohistochemical staining of human tissue showing IRF4 expression.