


Product Datasheet

EZH2 antibody (orb353501)

Description	EZH2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Histone-lysine N-methyltransferase EZH2 protein (1-250AA)
Target	EZH2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15910
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 (lane 1), J...