



Product Datasheet

SMARCC2 antibody (orb353503)


Description	SMARCC2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IHC, WB
Immunogen	Recombinant Human SWI/SNF complex subunit SMARCC2 protein (300-650AA)
Target	SMARCC2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TAQ2
Dilution Range	WB:1:200-1000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver tis...



Immunohistochemical staining of human to...



Immunoprecipitation analysis of Hela cel...