


Product Datasheet

ACTR1A antibody (orb353552)

Description	ACTR1A antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Alpha-centractin protein (187-376AA)
Target	ACTR1A
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P61163
Dilution Range	WB:1:500-5000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical
staining of human
lu...