



Product Datasheet

MRGPRX2 antibody (orb35358)

Description	MRGPRX2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Mas-related G-protein coupled receptor member X2 protein (286-330AA)
Target	MRGPRX2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96LB1
Dilution Range	WB:1:1000-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...



Immunofluorescence analysis of MCF-7 cel...