



Product Datasheet

IL18RAP antibody (orb353668)


Description	IL18RAP antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Interleukin-18 receptor accessory protein (390-599AA)
Target	IL18RAP
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95256
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human epidermis



Immunohistochemical staining of human lung tissue



Immunohistochemical staining of human colon tissue