


Product Datasheet

MAPKAP1 antibody (orb353714)

Description	MAPKAP1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Target of rapamycin complex 2 subunit MAPKAP1 protein (1-290AA)
Target	MAPKAP1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BPZ7
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa(Lane 1), M...