


Product Datasheet

ABLIM3 antibody (orb353763)

Description	ABLIM3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Actin-binding LIM protein 3 protein (266-544AA)
Target	Serine protease inhibitor Kazal-type 2-like protein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O94929
Dilution Range	WB:1:500-5000, IHC-P:1:200-500, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human adrenal gland tissue.