

Product Datasheet

ACADL antibody (FITC) (orb353773)

Description	ACADL antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Long-chain specific acyl-CoA dehydrogenase, mitochondrial protein (30-260AA)
Target	Pectate lyase 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P28330
Dilution Range	WB:1:500-5000, ICC:1:20-200, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.