

Product Datasheet

ACBD7 antibody (orb353775)

Description	ACBD7 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Acyl-CoA-binding domain-containing protein 7 protein (1-88AA)
Target	ACBD7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N6N7
Expiration Date	12 months from date of receipt.