

Product Datasheet

ACTR8 antibody (orb353799)

Description	ACTR8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Actin-related protein 8 protein (35-175AA)
Target	ACTR8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H981
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549(Lane 1), H...