



Product Datasheet

ACTRT1 antibody (orb353803)

Description	ACTRT1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Actin-related protein T1 protein (201-376AA)
Target	Pollen allergen Dac g 3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TDG2
Dilution Range	WB:1:500-5000, IHC-P:1:200-500, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human resected tissue showing ACTRT1 expression.



Immunohistochemical staining of human gastric tissue showing ACTRT1 expression.