

Product Datasheet

ACVR1B antibody (Biotin) (orb353816)

Description	ACVR1B antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Activin receptor type-1B protein (150-300AA)
Target	Outer capsid protein VP4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P36896
Dilution Range	WB:1:500-5000, IF:1:1-10, IP:1:200-2000
Expiration Date	12 months from date of receipt.