

Product Datasheet

ADAP2 antibody (orb353827)

Description	ADAP2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Arf-GAP with dual PH domain-containing protein 2 protein (110-381AA)
Target	ADAP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NPF8
Expiration Date	12 months from date of receipt.