

Product Datasheet

ANKRD7 antibody (orb353887)

Description	ANKRD7 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Ankyrin repeat domain-containing protein 7 protein (54-254AA)
Target	ANKRD7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q92527
Expiration Date	12 months from date of receipt.