

Product Datasheet

APOBEC3A antibody (Biotin) (orb353916)

Description	APOBEC3A antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human DNA dC->dU-editing enzyme APOBEC-3A protein (1-199AA)
Target	Beta-lactamase inhibitory protein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31941
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.