


Product Datasheet

Arg1 antibody (orb353923)

Description	Arg1 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IP
Immunogen	Recombinant Mouse Arginase-1 protein (1-323AA)
Target	Arg1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q61176
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of Rat live...