

Product Datasheet

Arg2 antibody (HRP) (orb353930)

Description	Arg2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Mouse Arginase-2, mitochondrial protein (23-354AA)
Target	Vertebrate ancient opsin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	008691
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.