

Product Datasheet

ARHGEF39 antibody (Biotin) (orb353940)

Description	ARHGEF39 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Rho guanine nucleotide exchange factor 39 protein (22-197AA)
Target	ARHGEF39
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N4T4
Expiration Date	12 months from date of receipt.