

Product Datasheet

ARMC8 antibody (Biotin) (orb353964)

Description	ARMC8 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Armadillo repeat-containing protein 8 protein (1-240AA)
Target	ARMC8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IUR7
Expiration Date	12 months from date of receipt.