

Product Datasheet

CISD2 antibody (HRP) (orb35397)

Description	CISD2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human CDGSH iron-sulfur domain-containing protein 2 protein (61-135AA)
Target	CISD2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N5K1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.