


Product Datasheet

ARRDC3 antibody (orb353983)

Description	ARRDC3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Arrestin domain-containing protein 3 protein (1-270AA)
Target	ARRDC3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96B67
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 whole cell...