



Product Datasheet

BRINP2 antibody (orb354071)

Description	BRINP2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human BMP/retinoic acid-inducible neural-specific protein 2 protein (34-300AA)
Target	BRINP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Antigen Affinity Purified
Uniprot ID	Q9C0B6
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sm...



Immunohistochemical staining of human co...