

Product Datasheet

C1orf123 antibody (Biotin) (orb354148)

Description	C1orf123 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human UPF0587 protein C1orf123 protein (1-160AA)
Target	C1orf123
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NWW4
Expiration Date	12 months from date of receipt.