

Product Datasheet

C5orf22 antibody (FITC) (orb354189)

Description	C5orf22 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human UPF0489 protein C5orf22 protein (101-420AA)
Target	C5orf22
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q49AR2
Expiration Date	12 months from date of receipt.