

Product Datasheet

Camk2n1 antibody (FITC) (orb354233)

Description	Camk2n1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Immunogen	Recombinant Mouse Calcium/calmodulin-dependent protein kinase II inhibitor 1 protein (1-78AA)
Target	Camk2n1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6QWF9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.