

Product Datasheet

CCDC181 antibody (FITC) (orb354285)

Description	CCDC181 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Coiled-coil domain-containing protein 181 protein (91-280AA)
Target	CCDC181
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5TID7
Expiration Date	12 months from date of receipt.