

Product Datasheet

CMTR2 antibody (Biotin) (orb354416)

Description	CMTR2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Cap-specific mRNA (nucleoside-2\'-O-)-methyltransferase 2 protein (301-600AA)
Target	CMTR2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IYT2
Expiration Date	12 months from date of receipt.