

Product Datasheet

COMMD9 antibody (FITC) (orb354445)

Description	COMMD9 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human COMM domain-containing protein 9 protein (1-198AA)
Target	COMMD9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9P000
Expiration Date	12 months from date of receipt.