

Product Datasheet

COX11 antibody (Biotin) (orb354452)

Description	COX11 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Cytochrome c oxidase assembly protein COX11, mitochondrial protein (115-276AA)
Target	COX11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y6N1
Expiration Date	12 months from date of receipt.