

Product Datasheet

CRYZL1 antibody (orb354479)

Description	CRYZL1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Quinone oxidoreductase-like protein 1 protein (1-87AA)
Target	CRYZL1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95825
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.

