



Product Datasheet

CSNK1D antibody (orb354491)


Description	CSNK1D antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Casein kinase I isoform delta protein (152-398AA)
Target	CSNK1D
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P48730
Dilution Range	WB:1:1000-5000, IHC-P:1:1000-2000, IF:1:50-500, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Sk...



Immunohistochemical staining of human ly...



Immunohistochemical staining of human li...