



Product Datasheet

Alpha E-catenin antibody (orb354495)


Description	Alpha E-catenin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Catenin alpha-1 protein (1-288AA)
Target	CTNNA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35221
Dilution Range	WB:1:500-2000, IHC-P:1:500-1000, IF:1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of human sm...